

Thermal Innovations and Products

# CPAC/HeatPAC Adapter Plates


▶ [Datasheet](#)

## CPAC/HeatPAC Adapter Plates

INHECO offers high precision adapters for the single heated and cooled positions (CPAC Microplate, CPAC Ultraflat, HeatPAC), to allow a perfect fit for all kinds of different disposables. The adapter plates are optimized for an excellent heat transfer to the disposables and enhance the uniformity over all wells and the heat up or cool down time. The installation and / or the exchange of adapter plates can be performed very easily within minutes.

### Microplate Adapters

Article#	Description	Pos.	Type	Volume	Disposable	Manufacturer	Disposable #	Comments
3200202	Flat Bottom Plate							Base plate for 7900016
3200203	96 pos. PCR Plate	96	PCR			Standard PCR Plates		
3200204	384 pos. PCR Plate	384	PCR			Standard PCR Plates		
3200279	Round Bottom Plate	96	Micro	355 µL	Greiner		65020x, 650261	U-Bottom
3200339	Round Bottom Plate	384	Micro	35 µL	Corning		3673-3679	U-Bottom, Round Well
3200374	Reservoir MWP		Reservoir	90ml	Seahorse		S30023, S30035	Flat Bottom
3200398	Round Bottom Plate	96	Micro	0,5ml	Axygen		P-450-V-C	V-Bottom
3200401	Square Well Plate	384	Micro	0,24ml	Greiner		781280, -87, -89	V-Bottom
3200402	Round Bottom Plate	96	Micro	0,34ml	Greiner		651201, -207, -209	V-Bottom
7900016	Flat Bottom Plate		Micro			Std.SBS Flat Bottom Plates		3200202 with position. frame



3200204: 384 Pos. PCR Plate



3200279: Round Bottom Plate



7900016: Flat Bottom Plate


### Tubes

Article#	Description	Pos.	Type	Volume	Disposable	Manufacturer	Disposable #	Comments
3200327	Deep Well 96	96	Tube			Micronic		Micronic Tubes
7900024	Tube Rack Multiprobe	22	Tube	20+1,5ml	Eppendorf			6x20ml, 16x1,5ml
7900029	Tube Rack 24x1,5ml	24	Tube	1,5ml	Eppendorf		1,5ml Tubes	Eppendorf Tubes
7900040	4 Troughs	4	Trough					
7900041	Falcon Tubes	8	Tube	50ml		BD Biosciences	352070	
7900045	Tube Rack 24x2,0ml	24	Tube	2,0ml		Sarstedt	72693005	

## Deep Well Adapters

Article#	Description	Pos.	Type	Volume	Disposable	Manufacturer	Disposable #	Comments
3200239	Round Well Adapter	96	Deep Well	1,0ml		Ritter	43001-0010/0016	Round Well, Riplate
3200240	Square Well Adapter	96	Deep Well	2,0ml		Ritter	43001-0020/0200	Square Well, Riplate
3200280	Deep Well 96	1	Reservoir	240ml		Axygen	RES-SW96-HP	V-Bottom, Square Well
3200281	Deep Well 96	96	Deep Well	2,2ml		Abgene	AB-0932	V-Bottom, Square Well
3200283	Masterblock 96	96	Deep Well	1,0ml		Greiner	78020x/215/26x	U-Bottom, Square Well
3200298	Flat Bottom 384	1	Half Deep	170ml		Seahorse	S30042	Flat Bottom, 384 Reservoir
3200306	Deep Well 96	96	Deep Well	2,0ml		Whatman	7701-5200, -05	U-Bottom, Square Well
3200311	Half Deep Well 384	384	Half Deep	240µL		Greiner	781270, 781281	V-Bottom, Square Well
3200313	Deep Well 96	96	Deep Well	1,0ml		Corning	3958, 3959	U-Bottom, Round Well
3200317	Deep Well 96	96	Deep Well	1,3ml		Nunc	26025x,278743/52	U-Bottom, Round Well
3200321	12 Columns Reservoir	12	Deep Well	21ml		Seahorse	S30019	Pyramid-Bottom, Square Well
3200326	12 Channel Reservoir	12	Deep Well	22ml		VWR	732-0585, -86	Square Well
3200337	DWP 96 High Temp.	96	Deep Well	1,2ml		Abgene	AB-0564	U-Bottom, Round Well
3200360	Deep Well 96	96	Deep Well	2,2ml		Qiagen	19585	U-Bottom, Square Well
3200361	DWP 96 High Temp.	96	Deep Well	2ml		Axygen	P-DW-20-C	U-Bottom, Round Well
3200363	Reservoir 384	384	Reservoir	92ml		Axygen	RES-SW384-LP	V-Bottom, Low Profile
3200364	8 Column Reservoir	8	Reservoir	32ml		Seahorse	S30020	Pyramid-Bottom
3200370	Deep Well 48	48	Deep Well	5ml		Seahorse	S30004	Pyramid-Bottom
3200371	4 Column Reservoir	4	Reservoir	73ml		Seahorse	S30051	Pyramid-Bottom
3200376	Half Deep Well 96	96	Half Deep	0,8ml		Abgene	AB-0765, -0859	V-Bottom
3200388	Deep Well 96	96	Deep Well	1,2ml		Trefflab	96.8996.9.01	U-Bottom
3200399	Deep Well 96 HT	96	Deep Well	1,3ml		Nunc	26025x,278743/52	U-Bottom, Square Well
3201140	Masterblock 96	96	Half Deep	0,5ml		Greiner	786201, 786261	V-Bottom, Square Well

Customized adapters are available on request.



3200283: Deep Well Masterblock 96



3200239: Round Well Adapter 96

**INHECO**  
Industrial Heating & Cooling GmbH

Fraunhoferstrasse 11  
82152 Martinsried / Munich  
Germany

Tel.: +49 (0)89 899 593 - 101  
Fax: +49 (0)89 899 593 - 149

E-mail: [sales@inheco.com](mailto:sales@inheco.com)  
[www.inheco.com](http://www.inheco.com)

This Datasheet can also be downloaded from our website. For further information, please contact our sales department.

