

Thermal Innovations and Products

# Shaker Adapter Plates

► [Datasheet](#)

## Shaker Adapter Plates

INHECO offers high precision adapters for the Variomag Teleshake and Thermoshake products, to allow a perfect fit for all kinds of different disposables. The adapter plates are optimized for an excellent heat transfer to the disposables and enhance the temperature uniformity over all wells and the heat up or cool down time. The installation and / or the exchange of adapter plates can be performed very easily within minutes.

### Microplate Adapters

Article#	Description	Pos.	Type	Volume	Disposable Manufacturer	Disposable #	Comments
3200245	Round Bottom Plate	96	Micro		Greiner		U-Bottom
3200292	Micro Plate 96	96	Micro	0,36ml	Corning	3799	U-Bottom
3200299	Micro Plate 96	96	Micro	0,45ml	Corning	3343	V-Bottom
3200300	Micro Plate 96	96	Micro	0,36ml	Corning	3355	U-Bottom
3200357	Micro Plate 384	384	Micro	0,04ml	Corning	367x	U-Bottom, Square Well
3200359	Micro Plate 384	384	Micro	0,14ml	Greiner	781280/287/289	V-Bottom, Square Well
3200386	Micro Plate 96	96	Micro	0,3ml	Nunc	249935	V-Bottom, Square Well
3200387	Micro Plate 96 PCR	96	PCR low		Standard PCR Plates		Low profile adapter
3201085	Micro Plate 96 PCR	96	PCR		Standard PCR Plates		
3201086	Micro Plate 384 PCR	384	PCR		Standard PCR Plates		
3201087	Flat Bottom Plate				Std. SBS Flat Bottom Plates		



3200245: Round Bottom Plate



3201085: 96 Pos. PCR Plate



3201087: Flat Bottom Plate

### Tubes

Article#	Description	Pos.	Type	Volume	Disposable Manufacturer	Disposable #	Comments
3200267	Tube Rack 96	96	Tube	0,5ml	Micronic		Micronic Tubes
3200394	Tube Rack 96x1,2ml	96	Tube	1,2ml	VWR	1,2ml Tubes	VWR Sample Library Tubes
7900031	Tube Rack 24x1,5ml	24	Tube	1,5ml	Eppendorf	1,5ml Tubes	Eppendorf Tubes
7900039	Tube Rack 24x1,5ml	24	Tube	1,5ml	Eppendorf	1,5ml Tubes	Eppendorf Tubes, Cap Holder



7900031: Tube Rack 24x1,5ml




3200267: Micronic Tubes

## Deep Well Adapters

Article#	Description	Pos.	Type	Volume	Disposable	Manufacturer	Disposable #	Comments
3200285	Masterblock 96	96	Deep Well	2,0ml		Greiner	78027x, 780285	V-Bottom, Square Well
3200286	24 Uniplate	24	Deep Well	10ml		Whatman	7701-5102/6102	U-Bottom, Square Well
3200290	Deep Well 96	96	Deep Well	1,0ml		Corning	3958	U-Bottom, Round Well
3200309	Deep Well 96	96	Deep Well	2,0ml		Corning	3960, 3961	V-Bottom, Square Well
3200333	Deep Well 96	96	Deep Well	1,0ml		Ritter	43001-0010, -16	U-Bottom, Round Well, riplate
3200334	Half Deep Well 96	96	Half Deep	0,5ml		Corning	3957	V-Bottom, Round Well
3200338	Deep Well 96	96	Deep Well	1,3ml		Nunc	26025x, 278743	U-Bottom, Square Well
3200349	Deep Well 96	96	Deep Well	2,2ml		Qiagen	19585	U-Bottom, Square Well
3200353	Deep Well 96	96	Deep Well	2,2ml		Abgene	AB0661	U-Bottom, Round Well
3200358	Half Deep Well 384	384	Half Deep	0,24ml		Nunc	269390	V-Bottom, Square Well
3200382	Deep Well 6x8	48	Deep Well	5,0ml		Seahorse	S30004	Pyramid Bottom
3200385	Deep Well 96	96	Deep Well	1,1ml		SPEware	SP1196	U-Bottom, Round Well
3200389	Deep Well 96	96	Deep Well	1,2ml		Trefflab	96.8996.9.01	U-Bottom
3200390	Deep Well 96	96	Deep Well	2,0ml		Axygen	P-DW-20-C	U-Bottom
3200391	RED Dialyse Plate	48	Deep Well			ThermoFisher	90085	U-Bottom, Flat Bottom
3200392	Deep Well 96	96	Deep Well	1,2ml		Abgene	AB-0564	U-Bottom, Round Well

Customized Adapters are available on request.




3200309: Deep Well 96



3200333: Deep Well 96



3200286: Uniplate 24



3200334: Half Deep Well 96

**INHECO**  
Industrial Heating & Cooling GmbH

Fraunhoferstrasse 11  
82152 Martinsried  
Germany

Tel.: +49 (0)89 899593-101  
Fax: +49 (0)89 899593-149

E-mail: [sales@inheco.com](mailto:sales@inheco.com)  
[www.inheco.com](http://www.inheco.com)

This Datasheet can also be downloaded from our website. For further information, please contact our sales department.

