



America

CERTIFICATE

No. U10 046515 0050 Rev. 01

Holder of Certificate: **INHECO**
Industrial Heating & Cooling GmbH
Fraunhoferstr. 11
82152 Martinsried
GERMANY

Certification Mark:



Product: **Laboratory Equipment**
(Thermal Cycler)

Tested according to:
UL 61010-1:2012/R:2024-11
UL 61010-2-010:2019/R:2023-09
UL 61010-2-101:2019/R:2023-12
CSA C22.2 No. 61010-1:2012/U4:2024-11
CSA C22.2 No. 61010-2-010:2019
CSA C22.2 No. 61010-2-101:2019

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. The certificate holder shall not transfer this certificate to third parties. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing, Certification, Validation and Verification Regulations (TCVVR)". For Canadian standards TÜV SÜD America Inc. is accredited by the Standards Council of Canada to ISO/IEC 17065.

Test report no.: 713300673

Date, 2025-09-08

(Ralph Fischer)

CERTIFICATE

No. U10 046515 0050 Rev. 01

Model(s):**Multi Tec Control (MTC) devices:**

Thermoshake
Thermoshake AC 2mm
Thermoshake AC 3mm
Thermoshake AC LC
Thermoshake AC 180
Teleshake AC
Teleshake AC USB 2mm
Teleshake AC USB 3mm
Teleshake 95
Teleshake 95 AC
Teleshake 95 AC 3mm
Teleshake 95 AC USB 2mm
Teleshake 95 AC USB 3mm
CPAC Microplate
CPAC Microplate HT 2-TEC
CPAC Ultraflat
CPAC Ultraflat HT 2-TEC
CPAC Ultraflat HT 2-TEC Down Vent
CPLC Ultraflat MTC
CPLC Ultraflat HT 2-TEC
HeatPAC
Heated Lid

Accessory for use with Multi Tec Control (MTC) devices:

Heat Exchanger Liquid cooling small

Brand Name(s):

inheco

inheco 



America

CERTIFICATE

No. U10 046515 0050 Rev. 01

Parameters:

Model	Rating
Thermoshake	24 Vdc / 5.3 A
Thermoshake AC	24 Vdc / 5.3 A
Thermoshake AC LC	24 Vdc / 5.2 A
Thermoshake AC 180	24 Vdc / 5.3 A
Teleshake AC 2mm	24 Vdc / 1.5 A
Teleshake AC 3mm	24 Vdc / 1.5 A
Teleshake AC USB 2mm	24 Vdc / 1.9 A
Teleshake AC USB 3mm	24 Vdc / 1.9 A
Teleshake 95	24 Vdc / 4.8 A
Teleshake 95 AC	24 Vdc / 5.3 A
Teleshake 95 AC 3mm	24 Vdc / 5.3 A
Teleshake 95 AC USB 2mm	24 Vdc / 6 A
Teleshake 95 AC USB 3mm	24 Vdc / 6 A
CPAC Microplate	12 Vdc / 5.0 A
CPAC Microplate HT 2-TEC	24 Vdc / 4.7 A
CPAC Ultraflat	12 Vdc / 4.5 A
CPAC Ultraflat HT 2-TEC	24 Vdc / 4.3 A
CPAC Ultraflat HT 2-TEC Down Vent	24 Vdc / 4.3 A
CPLC Ultraflat MTC	12 Vdc / 5.0 A
CPLC Ultraflat HT 2-TEC	24 Vdc / 4.08 A
HeatPAC	24 Vdc / 4.16 A
Heated Lid	24 Vdc / 4.8 A
Heat Exchanger Liquid cooling small	12Vdc / 50 W